

# Certificate of Conformity

|                                   |                                            |
|-----------------------------------|--------------------------------------------|
| <b>VWR Catalogue Number</b>       | 113-8380                                   |
| <b>Description</b>                | Chemical resistant overalls, white, medium |
| <b>Lot Number</b>                 | MT00367380                                 |
| <b>Country of Origin</b>          | Manufactured in China                      |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2022-07-21                                 |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

### Compliance to European Directives and Regulations:

|                                                                                                   |
|---------------------------------------------------------------------------------------------------|
| PPE Regulation (EU) 2016/425 on Personal Protective Equipment                                     |
| EC Type Examination Certificate N° GB18/961117 issued by Notified Body 0598, SGS Fimko Oy Finland |

### European Standards:

|                                                                                        |
|----------------------------------------------------------------------------------------|
| EN ISO 13982-1:2004 + A1:2010 - Protective clothing for use against solid particulates |
| EN 13034:2005 + A1:2009 - Protective clothing against liquid chemicals                 |
| EN 1149-5:2008 - Protective clothing - Electrostatic properties                        |
| EN 1073-2:2002 - Protective clothing against radioactive contamination                 |

**Material** Polyester yarn, rubber elastic, paper label

**Sterilization** Non-sterile.

|                                |                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ATP Assay</b>               | Not applicable.                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>DNase &amp; RNase Free</b>  | No information available.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>BSE/TSE</b>                 | The fabrics used to make the product may contain up to 500 ppm of tallow derived substances. Processing conditions of the raw materials used in manufacturing the fabrics meet or exceed the recommended conditions by Commission Regulation (EU) 722/2012, European Commission 2011/C 73/01 or ISO 22442-1 Annex C for protection against TSE (Transmissible Spongiform Encephalopathy)/BSE (Bovine Spongiform Encephalopathy). |
| <b>Cytotoxicity</b>            | Testing is conducted and have been shown to be non-toxic.                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Non-Pyrogenic Statement</b> | No information available.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Latex Statement</b>         | Both latex-free and latex-containing elastics are used. The latex-containing elastics in those garments and accessories are located in the waist elastic or in the thumb loops. The latex-containing elastics are covered by a stitching / covering thread to minimize the risk of direct skin contact with the elastic itself.                                                                                                  |
| <b>BPA Statement</b>           | Bisphenol are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                                                                                |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.                                                             |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.                                                                                                                                                                                                                  |
| <b>Storage Conditions</b>      | Store in ambient condition.                                                                                                                                                                                                                                                                                                                                                                                                      |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.