

Certificate of Conformity

VWR Catalogue Number	See List
Description	Reagent Bottle, Clear or Amber Soda-Lime Glass, Wide Neck with Black Screw Cap
Country of Origin	Manufactured in Germany
Date of Issue (yyyy-mm-dd)	2021-11-09

Clear Glass

Catalogue Number	Capacity [mL]	Thread	External Diameter [mm]	Height [mm]
215-4378	25	GL 28	36	57
215-4379	30	GL 32	36	68
215-4380	50	GL 32	44	75
215-4381	50	GL 40	44	69.5
215-4382	100	GL 40	50	92.5
215-4383	125	GL 40	54	98
215-4384	150	GL 45	57	104
215-4385	200	GL 55	63	109
215-4386	250	GL 55	69	113
215-4387	500	GL 55	84	154
215-4388	1000	GL 68	103.2	180.5

Amber Glass

Catalogue Number	Capacity [mL]	Thread	External Diameter [mm]	Height [mm]
215-4389	25	GL 28	36	57
215-4390	30	GL 32	36	68
215-4391	50	GL 32	44	75
215-4392	50	GL 40	44	69.5
215-4393	75	GL 40	48.5	79
215-4394	100	GL 40	50	92.5
215-4395	125	GL 40	54	98
215-4396	150	GL 45	57	104
215-4397	150	GL 55	63	109
215-4398	250	GL 55	69	113
215-4399	300	GL 55	70	127
215-4400	500	GL 55	84	154
215-4401	1000	GL 68	103.2	180.5
215-4402	2500	GL 80	138	281.5

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material

Clear or Amber Soda Lime Glass.
Screw Cap: Polypropylene (PP).

The material used to produce the products are listed in List of Authorised Substances of Regulation EU/10/2011.

The material grades conform to EU Directive 1935/2004 on materials and articles intended to come into contact with food.

Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	This product is free of cytotoxicity effects.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Diethylhexyl phthalate (DEHP) is not included in the manufacture or the formulation of the raw materials and are not expected to be present.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	No special storage condition required, avoid direct sunlight, dry (relative humidity <75 %).
Shelf Life	Not Applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature