

Data sheet

Article: Pasteur pipettes 150 mm, cotton plugged
Code no.: 925 02 01

Materials:

Glass

Soda-lime glass

Chemical properties:

- Hydrolytic resistance at 98°C: Class 3 acc. ISO 719
- Hydrolytic resistance at 121°C: Class 3 acc. ISO 720
- Acid resistance: acc. ISO 1776
- Alkali resistance: Class 2 acc. ISO 695

Plug

100% cotton bleached, chlorine-free

Standard:

ISO 7712

Hirschmann Laborgeräte
GmbH & Co.KG



Roland Benz
Quality Manager

Eberstadt, 15.02.2021