



Description: Sterile hypodermic needle for single use.

Components:

Description	Raw Materials
Hub	Polypropylene + master batch
Cannula	Stainless steel
Adhesive	Epoxy resin
Lubricant	Dimethyl silicone oil



Product Range:

Product Code	Colour	Cannula Ø		Cannula Length		Bevel	
AN*1838R1	Pink	18G	1.2 mm	1 1/2 inch	38 mm	Regular	11°
AN*1838S1	Pink	18G	1.2 mm	1 1/2 inch	38 mm	Short	17°
AN*1925R1	Ivory	19G	1.1 mm	1 inch	25 mm	Regular	11°
AN*1938R1	Ivory	19G	1.1 mm	1 1/2 inch	38 mm	Regular	11°
AN*1938S1	Ivory	19G	1.1 mm	1 1/2 inch	38 mm	Short	17°
AN*1950R1	Ivory	19G	1.1 mm	2 inch	50 mm	Regular	11°
AN*2025R1	Yellow	20G	0.9 mm	1 inch	25 mm	Regular	11°
AN*2038R1	Yellow	20G	0.9 mm	1 1/2 inch	38 mm	Regular	11°
AN*2116R1	Green	21G	0.8 mm	5/8 inch	16 mm	Regular	11°
AN*2125R1	Green	21G	0.8 mm	1 inch	25 mm	Regular	11°
AN*2138R1	Green	21G	0.8 mm	1 1/2 inch	38 mm	Regular	11°
AN*2150R1	Green	21G	0.8 mm	2 inch	50 mm	Regular	11°
AN*2232R1	Black	22G	0.7 mm	1 1/4 inch	32 mm	Regular	11°
AN*2238R1	Black	22G	0.7 mm	1 1/2 inch	38 mm	Regular	11°
AN*2316R1	Blue	23G	0.6 mm	5/8 inch	16 mm	Regular	11°
AN*2325R1	Blue	23G	0.6 mm	1 inch	25 mm	Regular	11°
AN*2332R1	Blue	23G	0.6 mm	1 1/4 inch	32 mm	Regular	11°
AN*2425R1	Purple	24G	0.55 mm	1 inch	25 mm	Regular	11°
AN*2516R1	Orange	25G	0.5 mm	5/8 inch	16 mm	Regular	11°
AN*2525R1	Orange	25G	0.5 mm	1 inch	25 mm	Regular	11°
AN*2613R1	Brown	26G	0.45 mm	1/2 inch	13 mm	Regular	11°
AN*2623R1	Brown	26G	0.45 mm	7/8 inch	23 mm	Regular	11°
AN*2716R1	Grey	27G	0.4 mm	5/8 inch	16 mm	Regular	11°
AN*2719R1	Grey	27G	0.4 mm	3/4 inch	19 mm	Regular	11°
AN*3013R1	Yellow	30G	0.3 mm	1/2 inch	13 mm	Regular	11°

Packaging:	Single unit:	Blister pack (film/paper)
	Unit box:	100 pieces Dimensions = 105 x 84 x 85 mm
		For item numbers AN*1950R1 and AN*2150R1: Dimensions = 105 x 84 x 110 mm
	Shipping carton:	10.000 pieces Dimensions = 450 x 440 x 450 mm
		For item numbers AN*1950R1 and AN*2150R1: 5.000 pieces Dimensions = 545 x 240 x 450 mm

Sterilization: With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Storage: Shelf life of 5 years
** The batch number and expiry date are printed on the packaging*
** Do not store at extreme temperature and humidity*

Precautions for Use: Check integrity of the unit package before use
 For single use only
 Do not use if unit package is damaged.
 Dispose of safely after single use to avoid risk of infection

Status of Regulation: Class IIa, sterile product following MDD 93/42/EEC

Notified body: TÜV Product Service GmbH, TÜV SUD Group 0123
 Manufacturer: Zhejiang Kindly Medical Devices Co. Ltd. (China)
 EC representative: Shanghai International Corp. GmbH (Europe) - Hamburg, Germany
 Distributor: TERUMO EUROPE N.V.
 Interleuvenlaan 40
 B - 3001 Leuven



Terumo Corporation
 44-1, 2-chome, Hatagaya
 Shibuya-ku, Tokyo 151-0072,
 Japan
 Tel.: +81 3 3374 8111
 Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
 Interleuvenlaan 40
 3001 Leuven, Belgium
 Tel.: +32 16 39 25 80
 Fax: +32 16 40 02 49

European Offices

**Terumo Europe N.V.
 Benelux Sales Division**
 Interleuvenlaan 40
 3001 Leuven, Belgium
 Tel.: +32 16 39 25 80
 Fax: +32 16 39 25 99
 The Netherlands:
 Tel.: 0800 0220396
 Fax: 0800 0220414

Laboratoires Terumo France S.A.
 Parc Ariane - Bât. Uranus
 1, rue Héliène Boucher
 78284 Guyancourt Cedex, France
 Tel.: +33 130 96 13 00
 Fax: +33 130 43 60 85

Terumo Deutschland GmbH
 Ludwig-Erhard-Strasse 6
 65760 Eschborn, Germany
 Tel.: +49 6196 80 230
 Fax: +49 6196 80 23 200

Terumo Deutschland GmbH
 Zweigniederlassung Spreitenbach
 Bodenackerstrasse 3
 8957 Spreitenbach, Switzerland
 Tel.: +41 56 419 10 10
 Fax: +41 56 419 10 11

Terumo Europe España SL
 Avda. Juan Carlos I, Nº 13 - 7ª Planta
 Edificio Torre La Garena
 28806 Alcalá de Henares (Madrid),
 Spain
 Tel.: +34 9021 01 298
 Fax: +34 9021 01 358

Terumo Italia SRL
 Via Simone Martini, 143/145
 00142 Roma RM,
 Italy
 Tel.: +39 0651 96 14 20
 Fax: +39 0650 30 407

Terumo Sweden AB
 Nya Varvet, Byggnad 90
 Sven Källfelts Gata 18
 426 71 Västra Frölunda,
 Sweden
 Tel.: +46 3174 85 880
 Fax: +46 3174 85 890

**Terumo Denmark
 Filial of Terumo Sweden AB**
 Tel.: +45 7020 93 80
 Fax: +45 7020 94 80

Terumo UK Ltd
 Tamesis
 The Causeway
 Egham, Surrey
 TW20 9AW, United Kingdom
 Tel.: +44 1784 476 220
 Fax: +44 1784 476 201

**Terumo Europe N.V.
 Moscow Representative Office**
 BC Northern Tower
 13Fr., 10 Testovskaya Street
 Moscow 123317, Russia
 Tel.: + 7 495 988 4740
 Fax: + 7 495 988 4739

**Terumo Corporation Dubai Branch
 (Middle East and Africa)**
 Al Masraf Tower, 22nd Floor
 P. O. Box 20291 Dubai
 U.A.E.
 Tel.: + 971 4 221 2220
 Fax: + 971 4 221 3330

**Addresses subject to change,
 please consult our Website.**

© Registered Trademark
 Published by Terumo Europe N.V.
 H0 131GB-0912TE-I(09.12)E