



VELTEK ASSOCIATES, INC.

PRODUCT RELEASE CERTIFICATE



ASEPTIC PROCESSING INC.

STERI-PEROX 6%

Date of Manufacture: 31-Jul-2025

Product STERI-PEROX 6%

Fill Size 16 oz

Cat No: SPER-16Z-6%-E

Lot #: 25-SPER-521566EC

Expiration: 01-Jul-2027

TEST	SPECIFICATION	RESULT
Heavy Metals (as Pb)	< 5 ppm	< 5 ppm
Assay	6.2 - 6.9%	6.7 %
Acidity	< 5.5 ml	2.4 mL
Non-Volatile Residue	< 30 mg	11 mg
Barium	No Turbidity	PASS
Preservative	< 50 mg	0.0 mg
LAL	< 0.25 EU/ml	PASS
Sterility		
USP Sterility	PASS	PASS

Ben Zilly

Authorized Signature - Quality Assurance

20 Aug 2025

Date

CONTAMINATION CONTROL INNOVATIONS SINCE 1981

Corporate Headquarters 15 Lee Boulevard, Malvern, PA 19355-1234 USA, T: 610-644-8335 F: 610-644-8336
www.sterile.com