

Certificate of Conformity

Catalogue Number	See list
Description	J.T.Baker Robotic tips, non-sterile, blistered
Country of Origin	Manufactured in Germany
Date of Issue (dd/mm/yyyy)	21/03/2022

Catalogue Number	Item Description	Volume
49000-0006	Conductive tips 300 µL tray	300 µL
49001-0106	Conductive tips 1100 µL filter tray	1100 µL
49002-0006	Conductive tips 200 µL S-tray	200 µL
49002-2006	Conductive tips 200 µL S-tray IVD	200 µL
49002-0106	Conductive tips 200 µL filter S-tray	200 µL
49002-2106	Conductive tips 200 µL filter S-tray IVD	200 µL
49002-0007	Conductive tips 200 µL C-tray	200 µL
49002-0107	Conductive tips 200 µL filter C-tray	200 µL
49003-0006	Conductive tips 1000 µL S-tray	1000 µL
49003-2006	Conductive tips 1000 µL S-tray IVD	1000 µL
49003-0106	Conductive tips 1000 µL filter S-tray	1000 µL
49003-2106	Conductive tips 1000 µL filter S-tray IVD	1000 µL
49003-0007	Conductive tips 1000 µL C-tray	1000 µL
49003-0107	Conductive tips 1000 µL filter C-tray	1000 µL
49033-0006	Conductive tips 1000 µL wide bore tray	1000 µL
49033-0106	Conductive tips 1000 µL wide bore filter tray	1000 µL
49006-0006	Conductive tips 50 µL S-tray	50 µL
49006-2006	Conductive tips 50 µL S-tray IVD	50 µL
49006-0106	Conductive tips 50 µL filter S-tray	50 µL
49006-2106	Conductive tips 50 µL filter S-tray IVD	50 µL
49006-0007	Conductive tips 50 µL C-tray	50 µL
49006-0107	Conductive tips 50 µL filter C-tray	50 µL
49008-0003	Conductive tips 300 µL tray	300 µL
49008-2003	Conductive tips 300 µL tray IVD	300 µL
49008-0105	Conductive tips 300 µL filter tray	300 µL
49008-2105	Conductive tips 300 µL filter tray IVD	300 µL
49009-0005	Conductive tips 1000 µL tray	1000 µL
49009-2005	Conductive tips 1000 µL tray IVD	1000 µL
49009-0105	Conductive tips 1000 µL filter tray	1000 µL
49009-2105	Conductive tips 1000 µL filter tray IVD	1000 µL
49010-0003	Conductive tips 50 µL tray	50 µL
49010-2003	Conductive tips 50 µL tray IVD	50 µL
49010-0105	Conductive tips 50 µL filter tray	50 µL
49010-2105	Conductive tips 50 µL filter tray IVD	50 µL
49018-0106	Conductive tips 50 µL QG filter tray	50 µL
49209-0105	Conductive tips MP 1000 µL filter tray	1000 µL
49209-2105	Conductive tips MP 1000 µL filter tray IVD	1000 µL
49102-0006	Non-conductive tips 200 µL tray	200 µL

49103-0006	Non-conductive tips 1000 µL tray	1000 µL
------------	----------------------------------	---------

Quality System Compliance

Products are manufactured in [ISO 13485](#) and [ISO 9001](#) certified plant. Products are inspected and controlled in accordance with the current applicable product specification.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Tips: Polypropylene Tray and S-tray: Polycarbonate or C-tray: Polypropylene Filter: Polyethylene
IVD	Robotic tips are available with EC conformity declaration as accessories for clinical diagnostic applications.
Manufacturing Environment	Robotic tips are manufactured under clean room conditions according to DIN EN ISO 14644-1.
Sterilization	Not applicable
ATP Assay	Tested with detection level $< 1 \cdot 10^{-12}$
PCR inhibitors Free	Product yield equivalent to control < 10 amplifiable targets
Human DNA Free	Tested with detection level < 1 pg
DNase & RNase Free	Tested with detection level $< 1 \cdot 10^{-6}$ (DNase) and $< 1 \cdot 10^{-9}$ (RNase)
Nucleic Acid Free	Not applicable
BSE/TSE	No use of any raw materials produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or
Cytotoxicity	Not applicable
Non-Pyrogenic Statement	Not applicable
Latex Statement	Not applicable
BPA Statement	Not applicable
DEHP Statement	Not applicable
RoHS	The raw materials used are compliant with the currently valid directive 2011/65/EU and 2002/95/EG (RoHS).

REACH Statement	The product does not contain any "SVHC" ingredients, in accordance with the currently valid regulation (EC) No 1907/2006 on REACH.
Material Statement	Not applicable
Storage Conditions	Keep dry and away from sunlight
Shelf Life	5 years

Disclaimer: Ritter states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.