



Product Specification Sheet

VWR SQUARE QUARTZ FLOW CELL, LIGHT PATH 10 mm

VWR Cat.No.	VWRI634-0695
Product Name	SQUARE.QUARTZ FLOW CELL (10mm)
Serial Number	634-0695
Material	ES Quartz Glass
Material Thickness	1.25±0.1 (mm)
Size	12.5*12.5*52 L*W*H (mm)
Optical Path Length	10±0.05 (mm)
Optical Path Width	4±0.1 (mm)
Cubic Measure	3.5 (ml)
Transmittance	≥79% (at 200nm)
Package	1pcs/box
Beam Height	15mm

VWR International bvba/sprl, Haasrode Research Park Zone 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium

<http://www.vwr.com>

Date of issue: Dec.03rd 2012