




VWR® Tube Stand

Description: Tube Stand, 7 Holes, Non-Sterile

Purpose: Suitable for 1.5ml, 2.0ml Micro Centrifuge Tube 15ml, 50ml Conical Centrifuge Tube

Materials: PP (Polypropylene) Color: Purple

Dimensions (Unit: mm)



Features

- 7 openings, suitable for holding 1.5mL and 2.0mL microcentrifuge tubes, 15mL and 50mL conical centrifuge tubes.
- DNase/RNase free and non-pyrogenic
- Autoclavable at 121°C and freezable to -80°C



VWR® Tube Stand

VWR EU EAN. No.	Description	Qty. per bag/case
525-1003	Tube Filter Unit Stand, 7 Holes, Suitable for 1.5ml, 2.0ml Micro Centrifuge Tube 15ml, 50ml Conical Centrifuge Tube, Non-Sterile	5/50