



Specification

1.03858.0500 EMB (Eosin-Methylene blue) agar acc. HOLT-HARRIS and TEAGUE
GranuCult® plus

Specification	
Appearance (clearness)	clear to slightly opalescent
Appearance (color)	reddish-brown to violet-brown
pH-value (25 °C)	6.9 - 7.3

Specification	
Growth (Escherichia coli ATCC 25922 (WDCM 00013))	good to very good
Growth (Escherichia coli ATCC 8739 (WDCM 00012))	good to very good
Growth (Salmonella typhimurium ATCC 14028 (WDCM 00031))	good to very good
Growth (Bacillus cereus ATCC 11778 (WDCM 00001))	none to poor
Metallic luster (Escherichia coli ATCC 25922 (WDCM 00013))	+
Metallic luster (Escherichia coli ATCC 8739 (WDCM 00012))	+
Metallic luster (Salmonella typhimurium ATCC 14028 (WDCM 00031))	-
Metallic luster (Bacillus cereus ATCC 11778 (WDCM 00001))	-

Incubation: 22 ± 2 hours at 35 ± 1 °C, aerobic

Stefanie Fischer
Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.