

Certificate of Conformity

VWR Catalogue Number	See List
Description	Racks for test tubes
Country of Origin	Manufactured in Pakistan
Date of Issue (yyyy-mm-dd)	2022-04-22

Polyamide Coated Steel

VWR Catalogue Number	Hole Dimension [mm]	No.of Holes	Width x Depth x Height [mm]
212-0713	18	12 (2 × 6)	130 × 50 × 70
212-0714	18	18 (3 × 6)	130 × 70 × 70
212-0715	18	24 (4 × 6)	130 × 90 × 70
212-0716	18	24 (2 × 12)	250 × 50 × 70
212-0717	18	36 (3 × 12)	250 × 70 × 70
212-0718	18	48 (4 × 12)	250 × 90 × 70
212-0719	18	100 (10 × 10)	210 × 210 × 70

Stainless Steel

VWR Catalogue Number	Hole Dimension [mm]	No.of Holes	Width x Depth x Height [mm]
212-0720	18	12 (2 × 6)	130 × 50 × 70
212-0721	18	18 (3 × 6)	130 × 70 × 70
212-0722	18	24 (4 × 6)	130 × 90 × 70
212-0723	18	24 (2 × 12)	250 × 50 × 70
212-0724	18	36 (3 × 12)	250 × 70 × 70
212-0725	18	48 (4 × 12)	250 × 90 × 70

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016 & F-899** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Stainless Steel or Polyamide coated Steel
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	Not applicable.
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above product is latex free.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Keep dry avoid Sunlight.
Shelf Life	5 Years.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.