

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Magnetic Wing Stirrers, Micro
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2026-05-28

Triangular

Catalogue Number	Height [mm]	Diameter [mm]	For	Pk
442-0495	16	10	Container with max. Ø 11 – 12 mm	10
442-0494	11	8	Container with max. Ø 9 – 10 mm	10
442-0493	9	5.5	Container with max. Ø 6 – 7 mm	10

Semicircular

Catalogue Number	Height [mm]	Diameter [mm]	For	Pk
442-0497	20	13	Container with max. Ø 14 – 15 mm	10
442-0496	10	13	Container with max. Ø 14 – 15 mm	10

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material

Virgin Polytetrafluoroethylene (PTFE).
Core: Alnico V permanent magnet.
Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI plastic (121 °C).

Temperature Resistance	-200 °C to +280 °C.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	DEHP are not used in the manufacture of the raw material and are not expected to be present.
RoHS	Product does not contain any of the restricted substance above the maximum concentration values as listed in Annex II of RoHS 2 Directive 2011/65/EU, as amended by EU 2015/863.
REACH Statement	The material used for the production of this container does not contain more than 0.1 % (w/w) of any substance of very high concern listed in the ECHA candidate list.
Storage Conditions	Ambient temperature.
Shelf Life	Shelf life is generally considered to be indefinite, within the confines of normal use and storage conditions.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

VWR International bv, Haasrode Researchpark 2020, Geldenaaksebaan 464, B-3001 Leuven, Belgium - <http://www.vwr.com>

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	04-2025
2	2.0	Modified USP Class VI information	05-2026