



Certificate of Conformity

VWR JOSEPH PAPER, GRADE 551

VWR International certifies that the products specified conform to all specifications and quality criteria mentioned in this document.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité mentionnés dans ce document.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle in diesem Dokument erwähnten Spezifikationen und Qualitätskriterien erfüllen.

Cat. No.	Description	L×W	Pk
111-5005	Plain	350 mm × 500 m	500
111-5006	Folded	350 mm × 500 m	500
111-5007	Plain	350 mm × 500 m	800
111-5009	Folded	350 mm × 500 m	40

Designation:

Joseph paper

Material: 100% cellulose nitrate

Binder: no



Certificate of Conformity

Characteristics	Values	Units	Standards
Weight	23	g/sqm	150536
Thickness	0.07	mm	
Filtration velocity		sec.	DIN 53137
		sec.	Herzberg
Pore size		µm	
Burst strength	0.2	kg/cm ²	interne
Ash content		%	DIN 53138
Capillary rise		mm	DIN 1508787
Capacity of absorption		µg/cm ²	
Permeability to air		l/min	100cm ² /200mmCE
Permeability to water		m ³ /m ² /h	
Tensile strength M.D.(Machine Direction)		N/5cm	
Tensile strength C.D.(Cross Direction)		N/5cm	

Applications:

Absorbent paper, suitable for cleaning and drying of glass containers, tubes, microscopic plates, bottles,

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **December 2017**