

OpenSciEd Biology Unit 3

REFILL KIT # 470357-856

Inheritance & Variation of Traits



Lesson # Included Items

8	p53 BRCA1 pedigree	1	pkg/35
9	Coupon - Wild type yeast, UV-sensitive yeast	1	
9	Toothpicks	60	
9	Transfer pipets, sterile, 1 mL	22	
9	Spreaders (replaces glass beads; see note below)	2	pkg/10
9	Yeast-extract dextrose medium	3	
9	Petri dish	1	pkg

Be careful to touch the spreader on the handle only. Using the short end of the spreader at an angle, open the petri dish slightly and spread the cell suspension thoroughly around to the edges of the plate.

WARD470357-794