

Certificate of Conformity

Product Code: LW4238

Item Description: Pipette Extended Fine Tip Sterile Pastette® Pack of 400 Single Wrap Use like Pasteur pipettes for liquid transfer Bulb draw 3.3ml Drop Volume 20µl

Lot Number: 22250765

Material: LDPE (virgin low-density polyethylene)

Sterility: Sterile

Method of Sterilisation: E-beam Irradiation

Where a Sterility Assurance Level is stated, it is ensured by irradiating the product according to ANSI/AAMI/ISO 11137-2006 guidelines.

Date of Manufacture: 27th June 2022

Expiry Date: 27th June 2026

Manufacturing Location: Mexico

Other information relating to this product:

- Length: 153mm
- Bulb Draw: 3.3ml
- Drops per ml: 50
- Drop Volume: 20µl
- No graduations

Alpha Laboratories Ltd hereby certifies that these products have been manufactured in a controlled environment in accordance with established manufacturing guidelines. Product inspection is achieved using appropriate AQL Levels per product approved requirements. All items have been inspected and meet packaging requisites.

Certified by: 
Robert Vint, Managing Director

Date: 02 May 2023