

Certificate of Conformity

VWR Catalogue Number	See List
Description	Micro Centrifuge Tube
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-04-28

Catalogue Number	Capacity [mL]	Colour	Sterility
525-1524	0.5	Natural	Non-Sterile
525-1525	0.5	Natural	Sterile
525-1527	0.5	Amber	Non-Sterile
525-1496	0.5	Amber	Sterile
525-1528	1.7	Natural	Non-Sterile
525-1529	1.7	Natural	Sterile
525-1497	1.7	-	Sterile
525-1531	1.7	Amber	Non-Sterile
525-1499	1.7	Amber	Sterile
525-1532	2.0	Natural	Non-Sterile
525-1533	2.0	Natural	Sterile
525-1535	2.0	Amber	Non-Sterile
525-1523	2.0	Amber	Sterile

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 13485:2016** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Material	Virgin Polypropylene (PP). Resins are USP Class-VI certified, FDA regulation CFR 21, and are FDA approved for food contact.
Sterilization	Product labelled as sterile are irradiated and dose released upon ISO 11137 recommended practices in effect at the time of validation meet a minimum requirement of 10 ⁻⁶ SAL (which is not applicable for Non-Sterile products).
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Endotoxin levels are tested by Limulus Amoebocyte Lysate (LAL) Assay. Product is free from Endotoxin.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	This product is free of Bisphenol A.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed.

REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Room Temperature.
Shelf Life	3 years from the date of manufacture (which is not applicable for Non-Sterile products).

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.