



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	961CS01
GENERAL DESCRIPTION	PHARMA REG.FLOCK.APP B/01PCS
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 50 film-film pouches each containing 1 regular flocked swab in tube, with red cap.

PRODUCT COMPONENTS

TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 155.00 mm ± 1.00 mm INSERTION DIAMETER: 9.00 mm ± 0.50 mm EXTERNAL TUBE DIAMETER: 12.00 mm ± 0.50 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm ± 0.20 mm CAP LENGTH: 32.80 mm ± 0.30 mm INSERTION DIAMETER: 9.30 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.04-11/12/2024) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	REGULAR FLOCKED SWAB	SPECIFICATIONS
		TOTAL LENGTH: 151.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm BREAKING POINT: 78.50 mm ± 1.00 mm FLOCKED LENGTH: 16.00 mm ± 2.00 mm FLOCKED DIAMETER: 5.50 mm ± 0.50 mm STICK MATERIAL: POLYSTYRENE FLOCKED MATERIAL: NYLON
CAP + TUBE	CAP + TUBE	SPECIFICATIONS
		TOTAL LENGTH: 175.50 mm ± 3.00 mm

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	500
PACKAGING	FILM-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	IONIZING IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	18 MONTHS

This bulletin (Rev.04-11/12/2024) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.