



Specification

1.10866.0500 WL nutrient agar for microbiology

Specification

Appearance (clearness)	clear
Appearance (color)	blue-green
pH-value (25 °C)	5.3 - 5.7

Specification

Growth (Candida albicans ATCC 10231 (WDCM 00054))	good to very good
Growth (Saccharomyces cerevisiae ATCC 9763 (WDCM 00058))	good to very good
Growth (Lactobacillus acidophilus ATCC 4356 (WDCM 00098))	fair to very good
Growth (Leuconostoc mesenteroides ATCC 9135 (WDCM 00108))	good to very good
Growth (Enterococcus faecalis ATCC 11700)	good to very good
Growth (Escherichia coli ATCC 11775 (WDCM 00090))	good to very good
Growth (Proteus mirabilis ATCC 29906 (WDCM 00023))	good to very good

Incubation: 4 days; 28 °C; aerobic.

Stefanie Fischer

Responsible laboratory manager quality control

This document has been produced electronically and is valid without a signature.