

Product Number: TX3290
Lot Number: 06100769
Manufactured Date: 3/5/2026
Manufacturing Location: USA

Certificate of Analysis


Product Description	
Test	Specification
Product Description	Sterile 70% Isopropanol, 1 gallon (3.8 l) bottle

Bacterial Endotoxins			
Test	Specification	Units	Result
Endotoxin	≤0.25 EU/mL	EU/mL	0.125

Chemical			
Test	Specification	Units	Result
Acidity	Not more than 1.0 mL of 0.020 N NaOH	mL	0.26
Appearance	Clear, colorless liquid	Not Applicable	Pass
Limit of Non-Volatile Residues	Does not exceed 5 mg (0.01%)	mg	0.10
Odor	Characteristic IPA	Not Applicable	Pass
Specific Gravity @ 20° C / 20° C	0.872 - 0.883	Not Applicable	0.881
Volume % v/v	Isopropyl Alcohol 68% - 72%	Percent	Pass

Expiration Date			
Test	Specification	Units	Result
Expiration Date	Month (#). Year (####)	Date	3.2028

Approved By:



Christopher D. Wall
Quality Assurance Manager

Date:

April 16, 2026

Time:

8:55:34 AM EST

TEX-QA-REC-8.2-001 Rev002

Confidential

Texwipe Americas
1210 South Park Drive
Kernersville, North Carolina 27284
United States
Tel: +1 336-996-7046
Toll free: +1 800-839-9473
Email: info@texwipe.com

Texwipe Europe/Middle East/Africa
ITW Contamination Control BV
Saffierlaan 5
2132 VZ Hoofddorp
The Netherlands
Tel: +31 88-1307-410
Email: info@itw-cc.com

Advanced Molding Co., Inc. | Texwipe Asia
No. 5 Circuit St., Light Industry & Science
Park 1,
Cabuyao City, Province of Laguna
Philippines 4025
Tel: +63 49-543-0241
Email: cs_asia@texwipe.com

ITW Contamination Control (Wujiang) Co., LTD
Texwipe China
No.4660 Pang Jin Road, Wujiang
Economic Development Zone,
Suzhou, Jiangsu, 215200 China
Tel: +86 512 6303 3700
Email: chinasales@texwipe.com