

Certificate of Conformity

VWR Catalogue Number	See List
Description	Disposable Weighing Boat
Country of Origin	Manufactured in Vietnam
Date of Issue (yyyy-mm-dd)	2026-02-10

Standard

Catalogue Number	Capacity [mL]	Dimension [mm]	Colour
611-0093	7	46 × 46 × 8	White
611-9168	7	46 × 46 × 8	Blue
611-9175	7	46 × 46 × 8	Black
611-0094	100	85 × 85 × 24	White
611-9169	100	85 × 85 × 24	Blue
611-9176	100	85 × 85 × 24	Black
611-0095	250	140 × 140 × 22	White
611-9170	250	140 × 140 × 22	Blue
611-9177	250	140 × 140 × 22	Black

Anti-static

Catalogue Number	Capacity [mL]	Dimension [mm]	Colour
611-9178	7	46 × 46 × 8	White
611-3113	7	46 × 46 × 8	Blue
611-9181	7	46 × 46 × 8	Black
611-9179	100	85 × 85 × 24	White
611-3114	100	85 × 85 × 24	Blue
611-9182	100	85 × 85 × 24	Black
611-9180	250	140 × 140 × 22	White
611-3115	250	140 × 140 × 22	Blue
611-9183	250	140 × 140 × 22	Black

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared that the products are inspected and controlled through whole-production processing in accordance with current applicable product specifications and quality control SOPs.

Inspection records are reviewed and signed off by qualified personnel for product release.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements, inspected in according to AQL Level II 4.0 sampling plan per MIL-STD-1916.

Product Specifications

Material	100% Virgin, Vacuum-molded of polystyrene homo-polymer with non-toxic pigment. Materials used meet FDA 21 CFR 177.1640 Compliance.
Resist temperatures	-10 to 70 °C.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.

RoHS	According to EU Directive 2015/863 No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Storage temperatures -20 °C to +60 °C (-4 °F to 140 °F).
Shelf Life	Indefinite Shelf Life.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	08-2021
2	1.0	Shelf Life Modified	01-2024
3	2.0	Modified Quality system compliance and disclaimer	07-2025
4	3.0	Country of Origin changed to Vietnam	02-2026