

Material No.: G193-04

:
Revision No.: 0

Product Specification

GMP Manufactured Product

Test	Specification
Appearance	Passes test
Refractive index	1,1 – 1,5 °
pH, For Information Only	
Viscosity at 25 °C, For Information Only	cP
Density	0,9 – 1,2 g/ml
Total Aerobic Microbial Count	≤ 100 CFU/ml
Total Yeasts and Molds	≤ 10 CFU/ml
Absorption maximum at 380–450 nm	

GMP Manufactured Product
For Laboratory, Research or Manufacturing use. Suitable for Use in Biotechnology

Packaging Site: Aurora Mfg Ctr & DC