


Recombinant E. coli RecA Protein (Functional) (A455)

Specifications:

Name:	Recombinant E. coli RecA Protein (Functional)
Applications:	SDS-PAGE, Functional Studies
Applications:	Studies on homologous recombination mechanism and SOS response. Useful in the screening using probe from library by promotion of DNA hybridization. Facilitate DNA observation by electron microscope due to nucleofilament formation with DNA.
Expression System:	Escherichia coli
Nature:	Recombinant
Protein Species:	Escherichia coli
Protein Length:	Full length protein
Molecular Weight:	38 kDa
Biologically Active:	Yes
Purity:	> 90% (by SDS-PAGE).
Product Form:	Liquid
Concentration:	1 mg/ml
Formulation:	Supplied in 20 mM Tris-HCl, pH 8.0, with 1 mM EDTA, 150 mM KCl, 1 mM DTT, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

Images:



Recombinant E. coli RecA Protein (Functional)