

Material	22555.460
Material description	Celite 512
Grade	GPR RECTAPUR
CAS Number	91053-39-3
Molecular formula	SiO ₂
Molecular mass	60.08

Characteristics	Specifications
pH (25°C; 10 %)	5.0 - 11.0
Wet density	300 - 473 g/l
Loss on ignition (1000°C)	≤ 0.4 %
Particle size (< 75 µm)	≥ 90 %

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Jessie Pijck, Head of Laboratory - Haasrode
VWR International bv; Geldenaaksebaan 464; BE-3001
Leuven; Belgium

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.