

Technical data sheet

Description

MicroTube Rack



Features

- For safe sample storage
- V bottom for optimal sample recovery
- Refillable alphanumeric marked rack

Plates


VWR Catalog Number	Description	Volume (mL)	Colour	Packed
45000-0000	Rack with 96 individual tubes	1.2	Clear/blue	10
45005-0000	Strips of 8 tubes	1.2	Clear/blue	10
45000-0008	Strip of 8 tubes	1.2	Clear	120
45003-0000	Individual tubes	1.2	Clear	960
45000-0006	Strip of 8 plug caps	-	Clear	120

Material

Microtube rack : PP

Technical drawing and measurements


VWR Catalog Number	Rack dimensions		
	Length (mm)	Width (mm)	Height (mm)
45000-0000	127.50	85.50	56.7
45005-0000	127.50	85.50	56.7
45000-0008	71.60	8.60	45.8
45003-0000	Ø 8.6	Ø 8.6	46.85
45000-0006	76.80	7.30	5.4



Microtube rack with 96 individual tubes



Strip of 8 plug caps



Strip of 8 tubes