

ORP buffer +475 mV, 250 ml

Specification Sheet (Part# 238322)

Product Specifications

Redox Value (mV)	475 mV
mV Accuracy	±5 mV
Stability (Months)	24
Certificate	No
Package	250 ml

© 2018 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

HAMILTON®

Web: www.hamiltoncompany.com

USA: +1-800-648-5950

Europe: +41 58 610 10 10

Hamilton Americas & Pacific Rim

Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 U.S.A.
Tel: +1-800-648-5950
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41 58 610 10 10
Fax: +41-58-610-00-10
contact@hamilton.ch

To find a representative in your area, please visit hamiltoncompany.com/contacts