



Product Specification Sheet

VWR COOLING POUCHES REFRIGERANT BRICKS

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	W×D×H	Weight	Pk
216-0188	114×102×38 mm	426 g	36
216-0189	102×229×38 mm	880 g	18
216-0190	235×286×25 mm	1590 g	12
216-0191	102×114×19 mm	199 g	72

Other Characteristics:

- Phase changes at 0 °C/32 °F
- It is a Semi-rigid foam.
- Returns to foam state after thawing.
- Keep their shape during all stages of use, reusable.
- It is suitable for use during storage and transportation.
- Each cooling pouch is vacuum-sealed for consistent weight and shape.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **September 2012**

