



## Certificate of Conformity

---

### VWR Standard Disposable Spatula

---

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

| Catalog No | Lot No   | Type     | Colour | Blade width (mm) | Length (mm) | Pk  |
|------------|----------|----------|--------|------------------|-------------|-----|
| 231-0103   | GD190019 | Standard | Blue   | 7                | 210         | 300 |

**Material:** 100% Virgin Food Grade Polypropylene

**BSE Statement:** No raw materials that contain, or are derived from, animals are specified in the manufacture of this product.

**Melamine Statement:** No Melamine is specified in the manufacture of this process.

**Latex Statement:** No Latex is specified in the manufacture of this process.

**FDA Statement:** The materials used, and the manufacturing process, meet FDA standards for plastics coming into contact with food.

Quality control for the product is performed upon receipt of each batch. Random samples of each batch are tested for adherence to size, weight and other manufacturing specifications; ability to withstand boiling, autoclaving and exposure to liquid nitrogen without damage to the product are also tested.

**Shelf Life:** If stored in a cool and dry location polypropylene resins should have an unlimited shelf life. Heat and direct sunlight will degrade most plastic polymers. Storage of the above referenced products out of direct sunlight and in normal lab conditions should result in an unlimited shelf life.

*This document has been produced electronically and is valid without a signature*

Date of issue: **December 2020**