





European Cat no	Product Description	Capacity (mL)	Brim Volume (ml)	Material (Bottle / Cap)	Bottle Shape	Sterile	Colour (Bottle / Cap)	Neck Size & Finish	Neck ID (mm)	Length (mm)	Width (mm)	Height Without Cap (mm)	Height With Cap (mm)	General Wall Thickness (mm)	Bottle Weight (Grm)	Cap Weight (Grm)	Gross Weight (Grm) ±0.5gm
215-4497	HDPE BOTTLE SQ 150 ml, Plug Cap RED-ST	150 mL	172	HDPE / PP	Square	Yes	Natural / RED	45 / 400	37.42	50.5	50.5	82.24	86.2	0.7mm	16	4.5	20.5
215-4498	HDPE BOTTLE SQ 250 ml, Plug Cap RED-ST	250 mL	272	HDPE / PP	Square	Yes	Natural / RED	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4499	HDPE BOTTLE SQ 500 ml, Plug Cap RED-ST	500ml	550	HDPE / PP	Square	Yes	Natural / RED	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4504	HDPE BOTTLE RECT 1000 ml, Plug Cap RED-ST	1000ml	1070	HDPE / PP	Rectangular	Yes	Natural / RED	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4505	HDPE BOTTLE RECT 2000 ml, Plug Cap RED-ST	2000ml	2120	HDPE / PP	Rectangular	Yes	Natural / RED	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4506	HDPE BOTTLE SQ 150 ml, EP liner Cap RED-ST	150 mL	172	HDPE / PP	Square	Yes	Natural / RED	45 / 400	37.42	50.5	50.5	82.24	86.2	0.7mm	16	4.5	20.5
215-4507	HDPE BOTTLE SQ 250 ml, EP liner Cap RED-ST	250 mL	272	HDPE / PP	Square	Yes	Natural / RED	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4508	HDPE BOTTLE SQ 500 ml, EP liner Cap RED-ST	500ml	550	HDPE / PP	Square	Yes	Natural / RED	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4509	HDPE BOTTLE RECT 1000 ml, EP liner Cap RED-ST	1000ml	1070	HDPE / PP	Rectangular	Yes	Natural / RED	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4510	HDPE BOTTLE RECT 2000 ml, EP liner Cap RED-ST	2000ml	2120	HDPE / PP	Rectangular	Yes	Natural / RED	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4511	HDPE BOTTLE SQ 150 ml, Plug Cap BLUE	150 mL	172	HDPE / PP	Square	NO	Natural / BLUE	45 / 400	37.42	50.5	50.5	82.24	86.2	0.7mm	16	4.5	20.5
215-4512	HDPE BOTTLE SQ 250 ml, Plug Cap BLUE	250 mL	272	HDPE / PP	Square	NO	Natural / BLUE	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4513	HDPE BOTTLE SQ 500 ml, Plug Cap BLUE	500ml	550	HDPE / PP	Square	NO	Natural / BLUE	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4514	HDPE BOTTLE RECT 1000 ml, Plug Cap BLUE	1000ml	1070	HDPE / PP	Rectangular	NO	Natural / BLUE	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4515	HDPE BOTTLE RECT 2000ml, Plug Cap BLUE	2000ml	2120	HDPE / PP	Rectangular	NO	Natural / BLUE	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4516	HDPE BOTTLE SQ 150 ml, EP liner Cap WHITE	150 mL	172	HDPE / PP	Square	NO	Natural / WHITE	45 / 400	37.42	50.5	50.5	82.24	86.2	0.7mm	16	4.5	20.5
215-4524	HDPE BOTTLE SQ 250 ml, EP liner Cap WHITE	250 mL	272	HDPE / PP	Square	NO	Natural / WHITE	45 / 400	37.42	57	57	100.24	102.2	0.7mm	20	4.5	24.5
215-4525	HDPE BOTTLE SQ 500 ml, EP liner Cap WHITE	500ml	550	HDPE / PP	Square	NO	Natural / WHITE	63 / 410	56.12	76.5	76.5	119.5	122.5	0.7mm	40	10.5	50.5
215-4526	HDPE BOTTLE RECT 1000 ml, EP liner Cap WHITE	1000ml	1070	HDPE / PP	Rectangular	NO	Natural / WHITE	63 / 410	56.12	94.9	76.3	178.2	181.18	0.7mm	55	10	65
215-4527	HDPE BOTTLE RECT 2000 ml, EP liner Cap WHITE	2000ml	2120	HDPE / PP	Rectangular	NO	Natural / WHITE	63 / 410	56.12	108	97	230	233	0.7mm	90	10	100
215-4650	HDPE BOTL RECT 2500 ML, CAP WHITE	2500ml	2700	HDPE / PP	Square	Yes	Natural / WHITE	83 / 2000	71.5	109	109	257	264.5	0.7mm	112	18	130
215-4651	HDPE BOTL RECT 2500 ML, CAP RED ST	2500ml	2700	HDPE / PP	Square	Yes	Natural / RED	83 / 2000	71.5	109	109	257	264.5	0.7mm	112	18	130
215-4652	HDPE AMB BOTL RECT 2700ML CAP YEL ST	2700ml	3100	HDPE / PP	Rectangular	Yes	Amber / Yellow	83 / 2000	71.5	125	115	250	253.36	0.7mm	122	18	140


Material

Bottle : HDPE, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.

Closure : Polypropylene, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed With Food grade colour.

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.



150 / 250 ml

500 ml

1 L / 2 L

2.5 L

2.7 L

Material Recycle Code



HDPE

LDPE

PP

Description: HDPE Bottle

Drawing:
Not to Scale



Release: 2025/04/03

Revision: 1