



Product Specification Sheet

VWR COVER GLASSES, PURE WHITE, HYDROLYTIC CLASS I

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Thickness	L×W	Pk
631-0119	0,08 mm	18×18 mm	1
631-0120	0,13 mm	18×18 mm	1
631-0121	0,16 mm	18×18 mm	1
631-0123	0,08 mm	22×22 mm	1
631-0124	0,13 mm	22×22 mm	1
631-0125	0,16 mm	22×22 mm	1
631-0126	0,19 mm	22×22 mm	1
631-0127	0,13 mm	24×24 mm	1

Raw material: D263 M Borosilicate glass

Compliant:

- Directive 98/79/CE on In-vitro-diagnostic medical devices
- Norm ISO 8255, part 1
- Norm DIN EN 980
- Norm ISO 15223

Storage:

- Keep dry
- Avoid large variations in temperature during both storage and usage. Cooling of the product may cause condensation and lead to condensed water forming between the glasses.
- Slide storage periods: see "use before" date

Other characteristics:

- Without bubbles or flaws
- Cleaned and polished
- Microscopic cover glasses for in-vitro-diagnostic use.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.