

Material	24360.368
Material description	Gelatine
Grade	Ph.Eur. powder
CAS Number	9000-70-8

Additional information

Characteristics	Specifications
Appearance	Conforms (see CoA/CoS)
Identification A	Passes test
Identification B	Passes test
pH (1 %)	3.8 - 7.6
Conductivity (30°C; 1 %; water)	≤ 1 mS/cm
Sulphur dioxide	≤ 50 ppm
Peroxides	≤ 10 ppm
Gel strength (Bloom value)	90 - 130
Fe (Iron)	≤ 30 ppm
Cr (Chromium)	≤ 10 ppm
Zn (Zinc)	≤ 30 ppm
Loss on drying (105°C)	≤ 15.0 %
TAMC	≤ 10 ³ CFU/g
TYMC	≤ 10 ² CFU/g
Escherichia coli	Absent
Salmonella	Absent
Residual solvents	Unlikely by manuf.process
Conforms to Ph.Eur.	Passes test

Signature	Additional information
We certify that this batch conforms to the specifications listed above. This document has been produced electronically and is valid without a signature. Jessie Pijck, Head of Laboratory - Haasrode VWR International bv; Geldenaaksebaan 464; BE-3001 Leuven; Belgium	Not for use as Active Pharmaceutical Ingredient (API)

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.