

# Certificate of Conformity

|                                      |                                                                                                                 |
|--------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>          | 221-0205                                                                                                        |
| <b>Description</b>                   | Funnels with Long Stem, 50 mm Funnel Outer Diameter, 10 mm Stem Outer Diameter, 85 mm Length, 45 mm Stem Length |
| <b>Lot Number</b>                    | J103076                                                                                                         |
| <b>Date of Manufacture (yyyy-mm)</b> | 2020-11                                                                                                         |
| <b>Country of Origin</b>             | Manufactured in United Kingdom                                                                                  |
| <b>Date of Issue (yyyy-mm-dd)</b>    | 2022-05-31                                                                                                      |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

## Product Specifications

|                               |                                                                                                                                                                                                                                                                                                  |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Virgin Polytetrafluoroethylene (PTFE) Grade TF1620.<br>Raw material is in compliance with Directive FDA 21 CFR 177.1550 and complies with USP class VI -27, NF 22,2004 for plastic at 70 °C.<br>The test is carried out by 3M-Dyneon who manufactures PTFE used and refers to PTFE raw material. |
| <b>Temperature Resistance</b> | -200 to +280 °C.                                                                                                                                                                                                                                                                                 |
| <b>Sterilization</b>          | Non-Sterile.                                                                                                                                                                                                                                                                                     |
| <b>ATP Assay</b>              | Not Applicable.                                                                                                                                                                                                                                                                                  |
| <b>DNase &amp; RNase Free</b> | Not Applicable.                                                                                                                                                                                                                                                                                  |

|                                |                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).          |
| <b>Cytotoxicity</b>            | Not Applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Non-Pyrogenic Statement</b> | Not Applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                                  |
| <b>BPA Statement</b>           | Bisphenol are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                    |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                   |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed. |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.                                                                                                                                                      |
| <b>Storage Conditions</b>      | No special storage conditions required. Generally clean environment, ambient temperature, free from the possibility of mechanical damage.                                                                                                                                                                                                                            |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.