

CERTIFICATE OF CONFORMANCE AND STERILITY

This certification is provided as full assurance that the following product code and lot number were manufactured in accordance via a Quality System certified to ISO 9001 with prescribed procedures and specifications.

Product Code:	MEQT0002SZ
Lot Number:	018543
Manufacturing Date:	2022.06.13
Expiration Date:	2024.06
Sterility Assurance Level:	10 ⁻⁶ (SAL determined by test method ISO11137)
Radiation Lot:	22062105
Irradiation Dose Specified:	25 kGy to 50 kGy
Minimum Absorbed Dose:	28.2 kGy
Maximum Absorbed Dose:	35.0 kGy
Radiation Date:	2022.06.21

Quality Department
Contec Cleanroom Technology (Suzhou) Co., Ltd.
2022-06-28



NORTH AMERICA

Contec, Inc.
525 Locust Grove
Spartanburg, SC 29303 USA
tel: +1 864 503 8333
fax: +1 864 503 8444

EUROPE

Contec Europe
ZI Du Prat, Avenue Paul Duplex
56000 VANNES France
tel: +33 (0)2 97 43 76 98
fax: +33 (0) 97 43 76 86

ASIA

Contec Cleanroom Technology (Suzhou) Co., Ltd.
17 Longyun Road
Suzhou Industrial Park, Suzhou 215024 China
tel: +86 512 6274 4050
fax: +86 512 6274 4051

web: www.contecinc.com
e-mail: info@contecinc.com

证书编号: Certificate No.: 0019739				辐照证明书 CERTIFICATE OF IRRADIATION	
执照登记号(Irradiation License No.): 91320582743739332H					
顾客: 康泰无尘科技(苏州)有限公司		编号: CHK/QR-8.2.6-02			
Customer:		版本/D/1			
产品名称 Product	包装规格(厘米) Dimension(cm)	产品批号 Product Lot No.	数量(箱) Qty.(ctns)	毛重(公斤) G.W. (kg)	
无尘拖把套MEQT0002SZ	——	018543	995 箱	——	
无尘包装袋BKT BAG	——	018542	244 箱	——	
涤纶针织无菌拖把头套MHQT0150	——	018555	200 箱	——	
无尘拖把头GCQT0150	——	018539	200 箱	——	
无尘包装袋BKT BAG2	——	018559	100 箱	——	
无尘拖把头TASK0250	——	018545	64 箱	——	
辐照日期 (Radiation 20220621)					
以下为空白					
张家港市中核华康辐照有限公司兹证明以上所列的所有产品(同制造商描述)在放射量测定系统精度和准确度的使用范围内接受到如下剂量。Zhangjiagang Municipal CNNC Huakang Radiation Co.,Ltd certifies that the products listed above cascades by the Manufacturer receiver the indicated doses within the precision and accuracy of the dosimetry system employed 剂量计测量系统的总不确定度为 5%					
辐照批号 Irradiation Lot No.	22062105	辐照完成日期 Irradiation finished at	20220622		
顾客要求最低吸收量 Specified minimum dose	25 kGy	实际最低吸收剂量 Minimum absorbed dose	28.2	kGy	
顾客要求最高吸收量 Specified maximum dose	50 kGy	实际最高吸收剂量 Maximum absorbed dose	35.0	kGy	
剂量计型号及批号 Dosimeter type & batch	210830	张家港市中核华康辐照有限公司 Zhangjiagang Municipal CNNC Huakang Radiation Co.,LTD 地址(ADD): 张家港市凤凰镇西张创业路28号 No.28,Chuangye Road,Fenghuang Town,Zhangjiagang City 电话(Tel) 0512-58421200 传真(Fax) 0512-58455118 邮编(P.C.) 215614			
检验员/日期: Inspected by:	何伟平 20220622				
批准人/日期: Approved by:	高爱民 20220622				
本证明书签字盖章有效		The certificate is valid after signature and stamped.			