

# Certificate of Conformity

|                                   |                                          |
|-----------------------------------|------------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                                 |
| <b>Description</b>                | VWR®, Inserts for Short Thread Vials ND9 |
| <b>Country of Origin</b>          | Manufactured in Germany                  |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2023-09-27                               |

| Catalogue Number | Description                                           | Capacity [mL] | Base           | Ø×H [mm] | Color |
|------------------|-------------------------------------------------------|---------------|----------------|----------|-------|
| 548-0812A        | Micro-insert with assembled plastic spring, silanized | 0.10          | Conical bottom | 5.7 x 29 | Clear |
| 548-3287A        | Micro-insert, silanized                               | 0.20          | Flat bottom    | 6 x 31   | Clear |

## Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

## Product Specifications

|                               |                                                                   |
|-------------------------------|-------------------------------------------------------------------|
| <b>Material</b>               | Clear glass, 1st hydrolytic class.                                |
| <b>Sterilization</b>          | Non-Sterile.                                                      |
| <b>ATP Assay</b>              | Not applicable.                                                   |
| <b>DNase &amp; RNase Free</b> | This product is free of any detectable DNase/RNase contamination. |

|                                |                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).                       |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Latex Statement</b>         | This product is not made from natural rubber latex.                                                                                                                                                                                                                                                                                                                                |
| <b>BPA Statement</b>           | Bisphenols are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                                 |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                                 |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr <sup>6+</sup> ), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analysis is performed. |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2023-06-14 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.                                                                                                                                                                    |
| <b>Storage Conditions</b>      | Dry, protected from light and standard warehouse temperature.                                                                                                                                                                                                                                                                                                                      |
| <b>Shelf Life</b>              | 2 years, conditions: store clean, dry, and protected from light at normal storage temperature.                                                                                                                                                                                                                                                                                     |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.

## Version History

| Sl.no | Version number | Description                          | Month & Year |
|-------|----------------|--------------------------------------|--------------|
| 1     | 1.0            | Created General CoC using new layout | 09-2023      |
|       |                |                                      |              |
|       |                |                                      |              |