

# Certificate of Conformity

|                                   |                                                                                       |
|-----------------------------------|---------------------------------------------------------------------------------------|
| <b>VWR Catalogue Number</b>       | See List                                                                              |
| <b>Description</b>                | Racks for microcentrifuge tube with 80 (5 × 16) holes,<br>225 × 67 × 29 mm Dimensions |
| <b>Country of Origin</b>          | Manufactured in China                                                                 |
| <b>Date of Issue (yyyy-mm-dd)</b> | 2021-11-24                                                                            |

| Catalogue Number | Colour           |
|------------------|------------------|
| 211-0204         | Natural          |
| 211-0205         | Orange           |
| 211-0206         | Yellow           |
| 211-0207         | Blue             |
| 211-0208         | Green            |
| 211-0209         | Pink             |
| 211-0210         | Assorted colours |

## Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

## QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements. inspected in according to AQL Level II 4.0 sampling plan per MIL-STD-1916.

## Product Specifications

|                               |                                                                                                                                         |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Material</b>               | Injection-molded of polypropylene homopolymer with pigments which are inert or minimally reactive with most common laboratory reagents. |
| <b>Sterilization</b>          | Not applicable.                                                                                                                         |
| <b>ATP Assay</b>              | Not applicable.                                                                                                                         |
| <b>DNase &amp; RNase Free</b> | Not applicable.                                                                                                                         |

|                                |                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>BSE/TSE</b>                 | No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).          |
| <b>Cytotoxicity</b>            | Not applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Non-Pyrogenic Statement</b> | Not applicable.                                                                                                                                                                                                                                                                                                                                                      |
| <b>Latex Statement</b>         | These products are not made from natural rubber latex.                                                                                                                                                                                                                                                                                                               |
| <b>BPA Statement</b>           | Bisphenol are not used in the manufacture of the raw material and are not expected to be present.                                                                                                                                                                                                                                                                    |
| <b>DEHP Statement</b>          | Please see RoHS information below.                                                                                                                                                                                                                                                                                                                                   |
| <b>RoHS</b>                    | No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed. |
| <b>REACH Statement</b>         | No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.                                                                                                                                                     |
| <b>Storage Conditions</b>      | Storage temperatures -20 °C to 60 °C (-4F to 140F).                                                                                                                                                                                                                                                                                                                  |
| <b>Shelf Life</b>              | Not applicable.                                                                                                                                                                                                                                                                                                                                                      |

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.  
This document has been produced electronically and is valid without a signature.