

Lab printers to fit your workflows

Need a new printer, but don't know which one? Brady has many printing solutions, each specifically designed to match your lab's volume, type of work, budget, connectivity and efficiency needs. Since sample identification is critical to lab work, it's essential to find the best fit. The following high-performance printers offer endless ways to boost your lab's efficiency – on demand and in step with your workflows.



**BMP®51
Label
Printer
(300 dpi)**

More mobility from the BMP®51 Label Printer lets you go where the work is: in the lab or in the field. With an intuitive QWERTY keyboard and backlit LCD display, this low-volume, handheld printer is your go-to for instant identification.

- ▶ Prints multi-line labels for the lab
- ▶ Large display and easy-to-use interface
- ▶ Ability to edit and print labels without being connected to a PC
- ▶ Ethernet (optional), USB, Wi-Fi (optional)
- ▶ Recommended materials: B-490, B-492, B-499



**BradyPrinter
M611 Mobile
Label Printer
(300 dpi)**

Design, edit and print sample labels right from your phone with the BradyPrinter M611 Mobile Label Printer. It's a mid-volume printer that uses the cloud to help you access files and import label data faster, and with more accuracy, anywhere you work.

- ▶ Prints using a mobile device via Bluetooth or Wi-Fi
- ▶ Allows label design and editing with Brady Workstation Express Labels Mobile App
- ▶ Built-in auto-cutter
- ▶ Bluetooth, USB 2.0, Wi-Fi
- ▶ Recommended materials: B-490, B-492, B-499



**BradyPrinter
i3300
Industrial
Label Printer
(300 dpi)**

Optimize lab workflows with greater speed using the i3300. It's a proven way to reduce label waste and sets up in a flash. This mid-volume industrial printer creates labels for vials, tubes, slides and all sample identification.

- ▶ Prints on the first label so there's no waste
- ▶ Foolproof material changeovers – labels and ribbons snap into place
- ▶ Handles a wide variety of lab identification
- ▶ Brady Workstation Laboratory ID Software Suite is available
- ▶ Ethernet, USB, USB 2.0, Wi-Fi
- ▶ Recommended materials: B-490, B-492, B-499



**BradyPrinter
i5300 Industrial
Label Printer
(300 or 600 dpi model)**

Get faster, more efficient precision printing for your lab with the i5300. Enjoy quick label setup and zero-waste material changeovers. The i5300 is an intuitive, high-mix, high-volume industrial printer for all lab applications.

- ▶ Prints small fonts with precision
- ▶ 20-second material changeovers with no-waste printing
- ▶ 7" color touchscreen and built-in auto-cutter
- ▶ Brady Workstation Laboratory ID Software Suite is available
- ▶ Cable, Ethernet, USB 2.0, USB Host
- ▶ Recommended materials: B-490, B-492, B-499

Need more options?

If you need another type of printer, software, barcode scanner or combo kit, reach out to see what's available. Remember, your lab is unique. Let's get the right fit for your workflow.



avantior™
delivered by **VWR™**

VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2022 Avantor, Inc. All rights reserved.

0322 Lit. No. 020146W