

Fire Assay Flux



VWR has more than 30 years of experience in manufacturing pre-mixed flux for fire assay fusion, achieving a global reputation for offering a homogeneous and high-quality product to the mining industry.

All our fluxes have low levels of gold content (typically <1ppb Au). Fluxes with low levels of gold and silver are also available (silver level is typically <0.17 ppm Ag). Both are available in fine powder or granular format with color varying between light yellow and orange.

Packaging is available in 20kg, 25kg, or 40kg UN-approved plastic pails or 500kg maxi-bags.

The product must be stored in a dry, indoor environment free from humidity.

VWR has manufactured custom flux recipes for decades and has the ability to satisfy the most demanding specifications.



VWR Mining US
Sparks, Nevada
Tel: 800-724-3620
Email: sparks@vwr.com

VWR Mining Canada
Mississauga, Ontario
Tel: 855-599-2692
Email: minescanada@vwr.com

VWR Mining Canada Exports
Vancouver, British Columbia
Tel: +1 (604) 270-2268
Email: vancouver@vwr.com