

Certificate of Conformity

VWR Pocket Thermometer

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Catalog No.	Accuracy (°C)	Range (°C)	Division (°C)	Length (mm)	Pk
620-0809	±1.0 °C	-5...+45	0.5	160	1
620-0810	±1.0 °C	-10...+110	1.0	160	1
620-0853	±1.0 °C	-35...+50	1.0	160	1

Material: Glass, Red Organic

Other Characteristics:

- These glass thermometers have a PTFE coating that helps to prevent contamination and contains liquid and glass in the event of breakage
- Red spirit. Not Traceable to NIST
- Total immersion thermometers.
- Ideal for field studies, wastewater treatment and any 'on the go' applications.
- Bulbous glass body.
- Individual serial number.
- Total Immersion, White Back Ground, Temperature range -35°C/50°C, Division 1
- Nylon window Case.

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **April 2019**