

Certificate of Conformity

VWR Catalogue Number	See List
Description	Storage boxes for sample vials and bottles, PP with transparent lid
Country of Origin	Manufactured in China - 548-1443A, 548-1444A, 548-1445A, 548-1446A, 548-1447A, 548-1448A Manufactured in Germany - 548-1995A, 548-1996A, 548-1997A, 548-1998A, 548-1999A, 548-2000A, 548-1511A, 548-1512A, 548-1513A
Date of Issue (yyyy-mm-dd)	2022-11-21

81 Cavities with Alphanumeric Coding

Catalogue Number	Colour	Dimension [mm]	Description
548-1443A	Blue	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials
548-1444A	Orange	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials
548-1445A	Pink	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials
548-1446A	Yellow	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials
548-1447A	Green	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials
548-1448A	Transparent	130 × 130 × 45	Box with 81 cavities for 1.5 ml vials

16 Cavities

Catalogue Number	Colour	Dimension [mm]	Description
548-1995A	Blue	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials
548-1996A	Orange	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials
548-1997A	Pink	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials
548-1998A	Yellow	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials
548-1999A	Green	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials

Catalogue Number	Colour	Dimension [mm]	Description
548-2000A	Transparent	67 × 67 × 45	Box with 16 cavities for 1.5 ml vials/2 ml shell vials

Other Boxes

Catalogue Number	Colour	Dimension [mm]	Description
548-1511A	Blue	130 × 130 × 102	Box with 25 cavities for 5/10/20 ml headspace vials
548-1512A	Violet	130 × 130 × 102	Box with 16 cavities for 20 ml EPA vial
548-1513A	Violet	130 × 130 × 102	Box with 10 cavities for 30/40 ml EPA vials

Quality System Compliance

Products are manufactured under the **ISO 9001:2015** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable product specifications.

Product Specifications

Material	Polypropylene.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from, or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.

Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
Temperature Resistance	-80 to 100 °C
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2022-06-10 are used in manufacturing of the raw materials used to produce this product. No routinely analyse is performed.
Storage Conditions	Not applicable.
Shelf Life	Not applicable.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.