




Europe Cat	Description	Capacity	Bottle Material	Cap Material	Vent Cap	Bottle Colour	Cap Colour	Cap Size	Outer Dia (A)	Height with Spout (B)	Height With Cap (C)	Spout Length (D)	Neck Inner Diameter (E)
VWR215-4294	Safety Wash Bottle Wide Mouth250 ml, Acetone	250ml	LDPE	PP	Yes	Natural	RED	38-415	Ø52	206	151	82	Ø27
VWR215-4295	Safety Wash Bottle Wide Mouth250 ml, Dichloromethane	250ml	LDPE	PP	Yes	Natural	Yellow	38-415	Ø52	206	151	82	Ø27
VWR215-4296	Safety Wash Bottle Wide Mouth250 ml, Distilled water	250ml	LDPE	PP	Yes	Natural	Natural	38-415	Ø52	206	151	82	Ø27
VWR215-4297	Safety Wash Bottle Wide Mouth250 ml, Ethanol	250ml	LDPE	PP	Yes	Natural	WHITE	38-415	Ø52	206	151	82	Ø27
VWR215-4298	Safety Wash Bottle Wide Mouth250 ml, Isopropanol	250ml	LDPE	PP	Yes	Natural	Yellow	38-415	Ø52	206	151	82	Ø27
VWR215-4299	Safety Wash Bottle Wide Mouth250 ml, Methanol	250ml	LDPE	PP	Yes	Natural	Green	38-415	Ø52	206	151	82	Ø27
VWR215-4300	Safety Wash Bottle Wide Mouth250 ml, Methyl ethyl ketone	250ml	LDPE	PP	Yes	Natural	Green	38-415	Ø52	206	151	82	Ø27
VWR215-4301	Safety Wash Bottle Wide Mouth250 ml, Sodium Hypochlorite (bleach)	250ml	LDPE	PP	Yes	WHITE	BLUE	38-415	Ø52	206	151	82	Ø27
VWR215-4302	Safety Wash Bottle Wide Mouth250 ml, Toluene	250ml	LDPE	PP	Yes	RED	RED	38-415	Ø52	206	151	82	Ø27
VWR215-4303	Safety Wash Bottle Wide Mouth250 ml, Write your Own Content	250ml	LDPE	PP	Yes	Natural	WHITE	38-415	Ø52	206	151	82	Ø27
VWR215-4304	Safety Wash Bottle Wide Mouth500 ml, Acetone	500ml	LDPE	PP	Yes	Natural	RED	38-415	Ø72	232	166	82	Ø44
VWR215-4305	Safety Wash Bottle Wide Mouth500 ml, Dichloromethane	500ml	LDPE	PP	Yes	Natural	Yellow	53-415	Ø72	232	166	82	Ø44
VWR215-4306	Safety Wash Bottle Wide Mouth500 ml, Distilled water	500ml	LDPE	PP	Yes	Natural	Natural	53-415	Ø72	232	166	82	Ø44
VWR215-4307	Safety Wash Bottle Wide Mouth500 ml, Ethanol	500ml	LDPE	PP	Yes	Natural	WHITE	53-415	Ø72	232	166	82	Ø44
VWR215-4308	Safety Wash Bottle Wide Mouth500 ml, Isopropanol	500ml	LDPE	PP	Yes	Natural	Yellow	53-415	Ø72	232	166	82	Ø44
VWR215-4309	Safety Wash Bottle Wide Mouth500 ml, Methanol	500ml	LDPE	PP	Yes	Natural	Green	53-415	Ø72	232	166	82	Ø44
VWR215-4310	Safety Wash Bottle Wide Mouth500 ml, Methyl ethyl ketone	500ml	LDPE	PP	Yes	Natural	Green	53-415	Ø72	232	166	82	Ø44
VWR215-4311	Safety Wash Bottle Wide Mouth500 ml, Sodium Hypochlorite (bleach)	500ml	LDPE	PP	Yes	WHITE	BLUE	53-415	Ø72	232	166	82	Ø44
VWR215-4312	Safety Wash Bottle Wide Mouth500 ml, Toluene	500ml	LDPE	PP	Yes	RED	RED	53-415	Ø72	232	166	82	Ø44
VWR215-4313	Safety Wash Bottle Wide Mouth500 ml, Write your Own Content	500ml	LDPE	PP	Yes	Natural	WHITE	53-415	Ø72	232	166	82	Ø44
VWR215-4314	Safety Wash Bottle Wide Mouth1000 ml, Acetone	1000ml	LDPE	PP	Yes	Natural	RED	53-415	Ø80	279	214	82	Ø44
VWR215-4315	Safety Wash Bottle Wide Mouth1000 ml, Dichloromethane	1000ml	LDPE	PP	Yes	Natural	Yellow	53-415	Ø80	279	214	82	Ø44
VWR215-4316	Safety Wash Bottle Wide Mouth1000 ml, Distilled water	1000ml	LDPE	PP	Yes	Natural	Natural	53-415	Ø80	279	214	82	Ø44
VWR215-4317	Safety Wash Bottle Wide Mouth1000 ml, Ethanol	1000ml	LDPE	PP	Yes	Natural	WHITE	53-415	Ø80	279	214	82	Ø44
VWR215-4318	Safety Wash Bottle Wide Mouth1000 ml, Isopropanol	1000ml	LDPE	PP	Yes	Natural	Yellow	53-415	Ø80	279	214	82	Ø44
VWR215-4319	Safety Wash Bottle Wide Mouth1000 ml, Methanol	1000ml	LDPE	PP	Yes	Natural	Green	53-415	Ø80	279	214	82	Ø44
VWR215-4320	Safety Wash Bottle Wide Mouth1000 ml, Methyl ethyl ketone	1000ml	LDPE	PP	Yes	Natural	Green	53-415	Ø80	279	214	82	Ø44
VWR215-4321	Safety Wash Bottle Wide Mouth1000 ml, Sodium Hypochlorite (bleach)	1000ml	LDPE	PP	Yes	WHITE	BLUE	53-415	Ø80	279	214	82	Ø44
VWR215-4322	Safety Wash Bottle Wide Mouth1000 ml, Toluene	1000ml	LDPE	PP	Yes	RED	RED	53-415	Ø80	279	214	82	Ø44
VWR215-4323	Safety Wash Bottle Wide Mouth1000 ml, Write your Own Content	1000ml	LDPE	PP	Yes	Natural	WHITE	53-415	Ø80	279	214	82	Ø44

Material : 100 % Virgin, FDA 21 CFR Compliance

BSE/TSE Statement: These product is manufactured without the use of Raw materials of animal origin.

Latex Free Statement: No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for this product.

Material Recycle Code



Description: Safety Wash Bottle

Drawing:
Not to Scale



Release: 6/17/2021

Revision: 0