

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Desiccators
Country of Origin	Manufactured in Czech Republic
Date of Issue (yyyy-mm-dd)	2026-01-07

Catalogue Number	Nominal size	Usable Height [mm]	Ext Diameter x H [mm]	Int Diameter [mm]	Pk
467-0069	100 DN	65	153 × 195	92	1
467-0070	150 DN	95	215 × 265	143	1
467-0071	200 DN	130	270 × 315	192	1
467-0072	250 DN	140	320 × 370	239	1
467-0073	300 DN	165	380 × 445	285	1

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations:

Material ISO 3585 Borosilicate glass 3.3, Clear.

Properties of Borosilicate Glass:

- Acid resistance Class I (to ISO 1776).
- Hydrolytic resistance Class I (HGB 1 to ISO 719, HGA1 to ISO 720).
- Alkali resistance Class II (to ISO 695).
- Coefficient of mean linear thermal expansion α : $3.3 \times 10^{-6} \text{ K}^{-1}$ (to ISO 7991; 20/300 °C).

Sterilization	Hot air sterilization, in the oven up to 140 °C. Steam sterilization, in an autoclave 121°C/ 20 min 134°C/ 10 min
ATP Assay	This product does not contain ATP.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	Product does not contain any substance from the REACH candidate list of SVHC.
Storage Conditions	Store at room temperature in normal warehouse conditions; Protect from any possible contamination.
Shelf Life	No limit.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	01-2026