

Product Specification

Meets Reagent Specifications for testing Ph. Eur. monographs

Test	Specification
Acidity or Alkalinity	passes test
Appearance of solution	passes test
Identification	passes test
Matter Insoluble in Water	passes test
Reducing Substances	passes test
Relative Density	0.790–0.793
Residue on Evaporation	≤ 50 ppm
Water (H ₂ O)	≤ 3 g/l
Additional Specification(s):	
Assay (by GC)	≥ 99.0 %
Related Substances:	
Benzene	≤ 2 ppm (v/v)
Methanol (CH ₃ OH)	≤ 0.05 % (v/v)
Other Impurities	passes test
2-Propanol (CH ₃ CHOHCH ₃)	≤ 0.05 % (v/v)

For Laboratory, Research or Manufacturing use.



Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000