

Biology

From the highest quality living and preserved specimens to the widest selection of media, microscope slides, and more, we have more options for your biology classroom than anywhere else.



Instruments and Tools Included in this Kit:

- 1 Teacher's Guide with MSDS
- 30 Student Worksheets and Guides
- 30 Microscopic slides
- 30 Cover slips
- 30 Tissue transfer brushes
- 30 Tissue treating trays with cavities A and B
- 1 Root tips scientifically grown and fixed, vial
- 4 Drop control bottles of Solution A (alcohol and hydrochloric acid)
- 4 Drop control bottles of Solution B (alcohol, acetic acid and chloroform)
- 4 Drop control bottles of Aceto Orcein Stain