

Hamilton Storage Solution

500 mL

Specification Sheet (Part/REF # 238931)



Hamilton storage solution is recommended to increase sensor longevity when not in use.

Product Specifications

Electrolyte	Hamilton Storage Solution, 500 mL
Package	500 mL
Stability Months (after date of manufacturing)	24
Storage Conditions	5 - 35 °C / 41 - 95 °F

Specifications are subject to change without notice
<https://www.hamiltoncompany.com/process-analytics/sensors/238931>

Spec. Version C