

Certificate of Conformity

VWR Catalogue Number	See List
Description	PP disposable beakers
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-07-07

VWR Catalogue Number	Capacity [mL]	Height [mm]	Ext. Diameter bottom (mm)	Ext. Diameter Top [mm]	Graduation Range [mL]
213-3915	15	38	28	37	5 - 15
213-3916	30	44	34	46	5-30
213-3917	50	52	39	52	5 - 50
213-3918	100	62	46	64	5 - 100
213-3919	150	79	49	70	5 - 150
213-3920	250	90	57	80	10 - 250
213-3921	400	111	68	96	25 - 400
213-3922	600	124	77	107	25 - 600

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material Polypropylene M800E.

Thickness 1 ± 0.12 mm.

Volume	600 mL ± 2.5 %
Sterilization	Non-Sterile.
PVC	Product is PVC free.
Phthalates	Product is phthalates free.
ATP Assay	Not applicable.
DNase & RNase Free	Not applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	Not applicable.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.
REACH Statement	Not applicable.
Storage Conditions	Room temperature.
Shelf Life	Not applicable.


Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Annex

Catalogue Number	Dimension [MM]								Volume
	A	B	C	D	E	F	R	T	
213-3922	106.10±0.58	123.40±0.64	1.0±0.12	77.50±0.46	1.50±0.12	14.00±0.18	8±0.16	1±0.12	600mL±2.5%
213-3921	94.70±0.52	110.50±0.58	1.0±0.12	68±0.46	1.50±0.12	13.90±0.18	6±0.16	1±0.12	400mL±2.5%
213-3920	80.40±0.52	90±0.52	1.0±0.12	56.60±0.46	1.50±0.12	13.90±0.18	3±0.14	1±0.12	250mL±2.5%
213-3919	70.20±0.46	79.80±0.46	0.80±0.12	49.30±0.36	1.20±0.12	13.60±0.18	3±0.14	0.80±0.12	150mL±2.5%
213-3918	63.30±0.40	62.30±0.46	0.65±0.12	45.80±0.36	1.20±0.12	10.50±0.18	3±0.14	0.80±0.12	100mL±2.5%
213-3917	53.20±0.40	52.70±0.40	0.65±0.12	38.70±0.32	1.20±0.12	12.60±0.18	2±0.12	0.65±0.12	50mL±2.5%
213-3916	45.30±0.36	44.40±0.36	0.65±0.12	33.70±0.32	1.20±0.12	12.60±0.18	2±0.12	0.65±0.12	30mL±2.5%
213-3915	37.40±0.32	38.50±0.36	0.65±0.12	28.20±0.28	1.20±0.12	12.60±0.18	2±0.12	0.65±0.12	15mL±2.5%



Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	02-2022
2	1.1	Added Annex	03-2022
3	1.2	Added PVC and Phthalates	06-2023
4	2.0	Modified Quality system compliance and disclaimer	07-2025