

Certificate of Conformity

VWR Catalogue Number	See List
Description	Cable ties, polyamide 6.6, VWR®
Country of Origin	Manufactured in Italy
Date of Issue (yyyy-mm-dd)	2022-05-13

Catalogue Number	Description	Colour	Dimension [mm]
129-1607	Cable ties, polyamide 6.6, resistance: 8	Black	100 × 2.5
129-1608	Cable ties, polyamide 6.6, resistance: 18	Black	150 × 3.6
129-1609	Cable ties, polyamide 6.6, resistance: 22	Black	200 × 4.8
129-1610	Cable ties, polyamide 6.6, resistance: 22	Black	300 × 4.8
129-1611	Cable ties, polyamide 6.6, resistance: 22	Black	370 × 4.8
129-1619	Cable ties, polyamide 6.6, resistance: 22	Blue	200 × 4.8
129-1620	Cable ties, polyamide 6.6, resistance: 22	Green	200 × 4.8
129-1621	Cable ties, polyamide 6.6, resistance: 22	Orange	200 × 4.8
129-1617	Cable ties, polyamide 6.6, resistance: 8	Red	100 × 2.5
129-1618	Cable ties, polyamide 6.6, resistance: 22	Red	200 × 4.8
129-1612	Cable ties, polyamide 6.6, resistance: 8	White	100 × 2.5
129-1613	Cable ties, polyamide 6.6, resistance: 18	White	150 × 3.6
129-1614	Cable ties, polyamide 6.6, resistance: 22	White	200 × 4.8
129-1615	Cable ties, polyamide 6.6, resistance: 22	White	300 × 4.8
129-1616	Cable ties, polyamide 6.6, resistance: 22	White	370 × 4.8
129-1622	Cable ties, polyamide 6.6, resistance: 22	Yellow	200 × 4.8

Quality System Compliance

Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Polyamide.
Operating temperature	–40 up to 85 °C.
Installation temperature	–10 up to 60 °C.
Melting point	260 °C.
Sterilization	Not applicable.
ATP Assay	Not applicable.
DNase & RNase Free	This product is free of any detectable DNase/RNase contamination.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not applicable.
Latex Statement	The above products are not made from natural rubber latex.
BPA Statement	Not applicable.
DEHP Statement	Not applicable.
RoHS	Not applicable.

REACH Statement No substances as listed on >Substance of Very High Concern< candidate list published 2022-01-17 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.

Storage Conditions No specific storage condition.

Shelf Life No shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature