

Certificate of Conformity

VWR Catalogue Number	See List
Description	Graduated Pipettes, Type 3, Class B, Borosilicate 9.1 glass
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2021-08-06

Catalogue Number	Capacity [mL]	Division [mL]	Tolerance [± ml]
612-4121	0.5	0.01	0.008
612-4122	1	0.01	0.01
612-4123	2	0.02	0.015
612-4124	5	0.1	0.045
612-4125	10	0.1	0.075
612-4126	20	0.1	0.15
612-4127	25	0.1	0.15
612-4128	50	0.2	0.3

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 17025:2017** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirement of the following standards and regulations:

Technical standards	DIN 12696
Calibration procedure	ISO 4787
Material	Borosilicate 9 Expansion glass as per ASTM E-438, type 2.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.