

Certificate of Conformity

VWR Catalogue Number	634-1212
Description	VWR® Absorption cuvettes with Teflon stopper and flat bottom, Quartz glass, Volume 700µl, Path length 2mm
Lot Number	2026041626
Date of Manufacture (yyyy-mm-dd)	2026-04-16
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2026-05-04

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with [GB/T19001-2016/ ISO 9001:2015](#).

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Path Length of Interior	:	2 +/- 0.03 mm
Path Width of Interior	:	10 +/- 0.20 mm
Path Length of Exterior	:	4.5 +/- 0.10 mm
Path Width of Exterior	:	12.5 +/- 0.10 mm
Height of Exterior	:	52 +/- 0.5 mm
Transmission (at 200nm)	:	> 80%

Product Specifications

Material	Optical glass or Quartz glass.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.

BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Not Applicable.
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Not Applicable.
Storage Conditions	Room temperature.
Shelf Life	Not Applicable.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacturer and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.