

**MP Biomedicals, LLC**

29525 Fountain Parkway
Solon, Ohio 44139

Telephone:
440/337 , 1200
Toll Free:
800/854 , 0530

Fax: 440/337 , 1180
web: www.mpbio.com

Certificate of Analysis

Product Description : Tryptose Phosphate Broth
Catalog Number : 1682149
Lot: C2487

Formula: N/A
CAS #: N/A
Physical Description : Clear, Light Yellow Solution

Formula Weight: N/A
Storage: 15 - 30° C

Test	Specification	Result
Sterility Test	Passes	Passes
pH	7.1 - 7.6	7.40
Osmolality	295 - 335 mM/kg	318 mM/kg
Functionality Test	Passes	Passes
Expiration Date	Month/Day/Year	10/31/2022

Mycoplasma: Negative

03/11/2020 - John Huang, PhD
MP Biomedicals, LLC.
Quality Control Manager

This is an electronically generated document
<mailto:biotech@mpbio.com>
<http://www.mpbio.com>

Online Ordering, MSDSs, certificates of analysis and data sheets now available on our web site
Technical Service: 1-800-279-5490 (440-337-1200) Customer Service: 1-800-854-0530 (440-337-1200)