

Certificate of Conformity

VWR Catalogue Number	See List
Description	Multi-layer adhesive floor mats, PureStep
Country of Origin	Manufactured in China
Date of Issue (yyyy-mm-dd)	2025-10-03

Adhesive film mats, PureStep, 30 layers, 8 µm ±20% mm thickness

Catalogue Number	L x W [mm]	Colour
115-0055	457 × 914	Blue
115-0056	457 × 914	Grey
115-0057	457 × 914	White
115-0061	457 × 1143	Blue
115-0062	457 × 1143	Grey
115-0063	457 × 1143	White
115-0047	610 × 762	Blue
115-0048	610 × 762	Clear
115-0049	610 × 762	Grey
115-0050	610 × 762	White
115-0067	610 × 914	Blue
115-0068	610 × 914	Grey
115-0069	610 × 914	White
115-0073	660 × 1143	Blue
115-0074	660 × 1143	Grey
115-0075	660 × 1143	White
115-0079	914 × 914	Blue
115-0080	914 × 914	Grey
115-0081	914 × 914	White
115-0085	914 × 1143	Blue

Catalogue Number	L x W [mm]	Colour
115-0086	914 × 1143	Grey
115-0087	914 × 1143	White
115-0091	914 × 1524	Blue
115-0092	914 × 1524	Grey
115-0093	914 × 1524	White

Adhesive film mats, PureStep, 60 layers, 8 µm ±20% mm thickness

Catalogue Number	LxW [mm]	Colour
115-0058	457 × 914	Blue
115-0059	457 × 914	Grey
115-0060	457 × 914	White
115-0064	457 × 1143	Blue
115-0065	457 × 1143	Grey
115-0066	457 × 1143	White
115-0051	610 × 762	Blue
115-0052	610 × 762	Clear
115-0053	610 × 762	Grey
115-0054	610 × 762	White
115-0070	610 × 914	Blue
115-0071	610 × 914	Grey
115-0072	610 × 914	White
115-0076	660 × 1143	Blue
115-0077	660 × 1143	Grey
115-0078	660 × 1143	White
115-0082	914 × 914	Blue
115-0083	914 × 914	Grey
115-0084	914 × 914	White
115-0088	914 × 1143	Blue
115-0089	914 × 1143	Grey
115-0090	914 × 1143	White

High tack mats, PureStep, 30 layers, 8 µm ±20% mm thickness

Catalogue Number	L x W [mm]	Colour
115-0167	457 × 1143	Blue
115-1774	457 × 1143	White
115-1775	660 × 1143	Blue
115-1776	660 × 1143	White
115-1748	914 × 914	Blue

High tack mats, PureStep, 60 layers, 8 µm ±20% mm thickness

Catalogue Number	L x W [mm]	Colour
115-0168	660 × 1143	White
115-0169	914 × 1143	White

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015 & ISO 14000:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Low Density Polyethylene, (LDPE) Blend, with 20% biohybrid resin.
Sterilization	Not Applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).

Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Tested.
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenols are not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	Please see RoHS information below.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Bis(2-Ethyhexyl) Phthalate (DEHP), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list published 2021-07-08 are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Store at 12-35 °C (54-95F).
Shelf Life	Maximum 5 years.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

Annex

Other Specifications:

Film : Low density Polyethylene

File Thickness/Layer : 1.8 mil / 50 micron

Adhesive : Acrylic

Adhesive Thickness : 0.3 mil/8 micron

Peel Adhesion : 6.5 oz./180 g/25mm +20 (individual layers)
24 oz./680 g/25mm +50 (bottom layers)

Tensile Strength: >1.2 kg/cm

Heat Resistance: 158° F/70° C -48 hours

Surface Tack Spec	Bottom Tack Spec
150-230	550-650

Other Characteristics:

- Anti-microbial, destroys germs and fungi
- Non-slip, colour-coded
- Reduced maintenance and cleaning costs

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	09-2021
2	1.1	Added Thickness	04-2022
3	1.2	Added other characteristic, tensile strength and heat resistance	09-2024
4	2.0	Modified Quality system compliance and disclaimer	07-2025
5	3.0	Modified shelf life minimum to maximum	10-2025