



Certificate of Conformity

VWR CLEANROOM PAPER, WITH SYNTHETIC LATEX, A3 BLUE 85 GSM, 5 REAMS OF 250 SHEETS, 1*1.250 ITEMS

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No. : 115-1952

Base Weight : 85 g/m²

Caliper : 97 ±3 µm

Tensile Strength : Machine Direction : 5.0 kg/15mm (Elmendorf Tear Test)
Cross Direction : 2.5 kg/15mm

Particle count (>0.5µm) : 170 pt/cm² (Liquid Particle Counter)

Anion :	Chloride (Cl)	0.420 µg/cm ²
	Nitrite (NO ₂)	<0.0015 µg/cm ²
	Bromide (Br)	<0.0015 µg/cm ²
	Phosphate (PO ₄)	<0.0038 µg/cm ²
	Sulfate (SO ₄)	<0.100 µg/cm ²

Cation :	Lithium (Li)	<0.0015 µg/cm
	Sodium (Na)	2.020 µg/cm
	Ammonium (NH ₄)	0.020 µg/cm
	Potassium (K)	0.100 µg/cm

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **February 2017**