



Certificate of Conformity

VWR Magnetic Stirring Bars, Boxed Sets

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht

Cat. No.	Type	Pk
442-0513	Octahedral 1	1 Set
442-0514	Oval 1	1 Set
442-0515	Polygon 1	1 Set
442-0516	Coloured Micro	1 Set
442-0517	Oval 2	1 Set
442-0518	Coloured Standard	1 Set
442-0519	Polygon 2	1 Set
442-0520	Mix 1	1 Set
442-0521	Mix 2	1 Set
442-0522	Cross-Shaped	1 Set

Core: Alnico V permanent magnet.

Compliant to: FDA and USP Class VI

PTFE used in the manufacture these products are pure, virgin material which conforms to FDA and USP Class VI requirements. It does not contain, or use in its manufacture any Class 1 ODC according to the Clean Air Act 1990, PFOA, POP, DEHIP and BPA, no animal material is used in the manufacture or is present in this material. PTFE used conforms to current REACH requirements.

Set Octahedral 1, PTFE covered, 12 octagonal magnetic stirring rods with centre ring, in practical, transparent box. Set includes two stirrers with each of following dimensions, L×Ø (mm): 13×8, 15×8, 25×10, 38×10, 51×10 and 64×10.

Set Cross-Shaped, PTFE covered, five cross-shaped magnetic stirring rods in practical, transparent box. Set includes one stirrer with each of following dimensions, L×Ø (mm): 10×10, 20×20, 25×25, 30×30 and 38×38.



Certificate of Conformity

Set Oval 1, PTFE covered, 12 oval magnetic stirring rods in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 10×5, 15×6, 25×10, 30×10 and 35×13.

Set also includes one stirrer with each of following dimensions, L×Ø (mm): 50×17 and 70×27.

Set Oval 2, PTFE covered, 13 oval magnetic stirring rods in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 20×10, 25×12, 30×16, 35×16, 40×20 and 50×20. Set also includes one stirrer with following dimensions, L×Ø (mm): 70×20.

Set Coloured Micro, PTFE covered, 12 coloured micro magnetic stirring rods in practical, transparent box.

Set includes one stirrer in each colour (blue/red/yellow) with following dimensions L×Ø (mm): 8×1.5, 5×2, 7×2 and 10×3.

Set Coloured Standard, PTFE covered, 24 coloured octagonal magnetic stirring rods with centre ring, in practical, transparent box.

Set includes one stirrer in each colour (blue/red/yellow) with following dimensions L×Ø (mm): 13×3, 13×8, 15×8, 22×8, 25×8, 38×8, 51×8 and 75×13.

Set Polygon 1, PTFE covered, 18 cylindrical magnetic stirring bars in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 10×6, 15×6, 20×7, 25×8, 30×8, 40×8, 50×8, 60×10 and 80×10.

Set Polygon 2, PTFE covered, 18 cylindrical magnetic stirring bars in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 10×6, 15×4.5, 20×6, 25×6, 30×6, 40×8, 50×8, 60×10 and 70×10.

Set Mix 1, PTFE covered, 22 cylindrical magnetic stirring bars, 16 with removable pivot ring, in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 13×8, 13×10, 20×8, 25×8, 25×10, 38×8, 42×10, 51×8 and 51×10. Set also includes one stirrer with each of following dimensions, L×Ø (mm): 32×8, 32×10, 38×10 and 45×8.

Set Mix 2, PTFE covered. Set of 14 magnetic stirring bars: 4 micro and 10 octagonal with centre ring, in practical, transparent box.

Set includes two stirrers with each of following dimensions, L×Ø (mm): 7×2, 10×3, 13×8, 15×8, 25×10, 38×10 and 51×10.

This document has been produced electronically and is valid without a signature

Disclaimer: VWR wants to state that this declaration will not discharge the user from his obligation to ensure the product is suitable for the intended use. - If there is any doubt ask us.

Date of issue: **September 2019**